

Administering Therapy Using a Freedom60® Syringe



1. Prepare a clean work area. Wash hands thoroughly.
2. Be sure medication is room temperature before you begin filling the syringe.
3. Verify you have the correct tubing and needle set.
4. Remove the sterile cap from the end of the tubing set and connect the luer disc end of this tubing to syringe.
5. Wind the large knob clockwise to move black tab to end of track. With syringe markings facing up load syringe into pump, making sure the luer disc is fully seated in end of pump.
6. Prepare sites, insert needles and cover with adhesive dressing. Begin infusion.

